



BREAS MEDICAL INTRODUCES THE VIVO 50 AND VIVO 60 “PEARL” VERSION

MOUThPIECE VENTILATION AND OTHER ENHANCEMENTS FOR IMPROVED TREATMENT OPTIONS AND PATIENT COMFORT IN HOMECARE VENTILATION THERAPY.

The highly successful Vivo 50 and Vivo 60 homecare ventilators are now further enhanced with a dedicated mouthpiece ventilation mode, upgraded performance and improvements of the Service Software. A complete package, called the “Pearl” version.

“We listen to patients and customers and continuously develop our products”, says Damian deGoo, CEO for Breas Medical. “We now introduce a completely new version of our highly successful Vivo 50 and Vivo 60 ventilators to meet demands of customers and to further enhance treatment options and patient comfort. We call it the ‘Pearl’ version.”

What are the benefits of the ‘Pearl version’?

The main feature of the ‘Pearl version’ is the new mouthpiece ventilation (MPV) mode. Developed by the Breas R&D team, this MPV mode offers some unique advantages such as no flow between breaths, a settable apnea alarm, an adapted low-pressure alarm and clear information messages on the display. “Our advanced algorithm solves some of the drawbacks with MPV treatment that clinicians informed us about”, says Roger Nyström, product manager at Breas. “Our MPV mode offers some clear advantages to patients and clinicians”. The performance upgrade consists of a settable breath rate of 0 in MPV mode and the display of Vte and MVe when using a leak circuit. The improvements of the service software help increase ease of use and save time for technicians. Altogether, the ‘Pearl’ version contains a set of enhancements that further widens the usage of the Vivo 50 and Vivo 60 ventilators and meet the growing demand for mouthpiece ventilation treatment.

When will the Vivo 50 and Vivo 60 ‘Pearl’ be available?

Vivo 50 and Vivo 60 ventilators with the ‘Pearl’ version will become available by the end of June. The release will be phased

and the introduction may be subject to local registration procedures. Vivo 50 and Vivo 60 ventilators in the field can be upgraded to the ‘Pearl’ version to make all new benefits also available to current patients and caregivers.

Breas Medical

Breas Medical is one of the global leaders in Respiratory and Sleep therapy. Breas’ products combine attractive Scandinavian design, excellent quality and user friendliness to enable patients with respiratory diseases to be transitioned quickly from the hospital to home. Breas Medical operates in more than 40 countries through a network of subsidiaries and distributors. Breas Medical employs approximately 150 people and successfully markets the Vivo, iSleep, Nippy, Clearway and Z1 brands. For more information, visit www.breas.com.

Vivo 50 and Vivo 60 homecare ventilators

The Vivo 50 and Vivo 60 ventilators are life-support ventilators for use in the hospital and home environment. The devices are recognized and appreciated for their combination of outstanding quality, intuitive user interface, extensive monitoring possibilities and robust design. The Vivo 60 is especially designed for use with pediatric in patients from 5 kg. <http://breas.com/products/>

Contact

Sebastian Mommers, VP Marketing,
sebastian.mommers@breas.com, +46 706 229908